



permobil

TILITE MANUAL WHEELCHAIRS

AERO Z

STANDARD IS NOT
IN OUR DNA.

This is where TiFit begins. The aluminium Aero Z frame is built perfectly around you. We utilise the most advanced computer-aided designs to craft a chair that optimises your fit and performance. TiFit... not standard.



TIFIT. PERSONAL ERGONOMICS

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

CUSTOM CONFIGURATION

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options than any other manufacturer!

TECHNICAL SPECIFICATIONS

Frame Features	Aluminium Mono-Tube Rigid TiFit Full Adjustability
Total Chair Width*	23" (58.4cm)
Seat Width	10" – 20" (25.4 – 50.8cm)
Seat Depth	10" – 20" (25.4 – 50.8cm)
Front Seat Height	16" – 21" (40.6 – 53.3cm)
Rear Seat Height	12.5" – 20" (31.7 – 50.8cm)
Front Frame Angle	70°, 75°, 80°, 85°, or 90°
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" (21.5 – 53.3cm) Adj. Angle 80° – 101°
Camber	0°, 2°, 4°, 6°, 8° or 12°

Center of Gravity Adjustment	6" (15.2cm)
Front Seat Width (Seat Taper)	0", 1.5", 2.5" (0, 3.8, 6.3cm)
Transit Options	Transit Tie-Down
Transport Weight*	5.6 kg
Weight Limit	120.2 kg
Warranty	Five Years on Frame One Year on TiLite Components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.2 kg.


COLOUR OPTIONS

HIGH GLOSS PAINT



TiLite Red Burgundy Neon Pink Super Wet Black

MATTE PAINT



Deep Blue Ice Titanium Grey Black Lava

CANDY PAINT



Ruby Red Azure Blue Electric Plum

PEARLESCENT PAINT



Midnight Blue Pearl White Pearl

METALLIC PAINT



Salsa Red Metallic Tangerine Metallic Acid Green Metallic Ocean Blue Metallic

ANODIZED COLOURS



Red Pepper Burnt Orange Sapphire Ultraviolet Fuchsia Black