



COMPARE AND CHOOSE

Your electric **handbike**


BATEC Mini 2

BATEC Hybrid

BATEC Electric 2

BATEC Rapid 2

BATEC Scrambler



BATEC MINI 2



BATEC HYBRID



BATEC ELECTRIC 2



BATEC RAPID 2



BATEC SCRAMBLER



Maximum speed	Maximum speed	Maximum speed	Maximum speed	Maximum speed
20 Km/h	N/A	25 Km/h	30Km/h	30 Km/h
Autonomy	Autonomy	Autonomy	Autonomy	Autonomy
28 km * with a 280 Wh battery 52 km * with a 522 Wh battery	27-40km * with a 280 Wh battery 50-75km * with a 522 Wh battery	22 km * with a 280 Wh battery 40 km * with a 522 Wh battery	22 km * with a 280 Wh battery 50 km * with a 678 Wh battery	22 km * with a 280 Wh battery 50 km * with a 678 Wh battery
Motor	Motor	Motor	Motor	Motor
350 RPM high-torque brushless motor Nominal power: 350 w Peak power: 600 w Reverse gear	200 RPM high-torque brushless motor Nominal power: 900 w	250 RPM high-torque brushless motor Nominal power: 900 W Peak power: 1050 W Reverse gear	300 RPM high-torque brushless motor Nominal power: 1200 W Peak power: 1300 W Reverse gear	300 RPM high-torque brushless motor Nominal power: 1200 W Peak power: 1300 W Reverse gear
Battery and charger	Battery and charger	Battery and charger	Battery and charger	Battery and charger
2 options: 36 V - 280 Wh or 36 V - 522 Wh 4A Charger	2 options: 36 V - 280 Wh or 36 V - 522 Wh 4A Charger	2 options: 36 V - 280W h or 36 V - 522 Wh 4A Charger	2 options: 36/48 V - 280W h or 48 V - 678 Wh 3A Charger	2 options: 36/48 V - 280W h or 48 V - 678 Wh 3A Charger
Frame	Frame	Frame	Frame	Frame
7005 T6 Aluminium/Aluminium 6061 T6 Black 1.1/8 headset Folding stem bar Removable right dropout Folding handlebar Matte black aluminium handlebar Available for quadriplegics and hemiplegics	7005 T6 Aluminium 1.1/8 FSA Orbit MX Black steering Stand frame with suspension - 170mm Aluminium Crank Arms - Available for quadriplegics and hemiplegics	7005 T6 Aluminium Black 1.1/8 headset - Removable right dropout - Matte black aluminium handlebar Available for quadriplegics and hemiplegics	7005 T6 Aluminium 1.1/8 FSA Orbit Mx Black steering - Removable right dropout - Matte black aluminium handlebar Available for quadriplegics and hemiplegics	7005 T6 Aluminium 1.1/8 FSA Orbit Mx Black steering - Removable right dropout Stand frame with suspension Matte black aluminium handlebar Available for quadriplegics and hemiplegics
Brakes	Brakes	Brakes	Brakes	Brakes
Avid BB5 mechanical brakes Brake levers with parking brake Avid G2 140-mm rotor disc	Avid BB7 mechanical brakes Avid speed dial 7 adjustable brake levers 8 Gear Shimano Alivio Gear Shifter	Avid BB5 mechanical brakes Brake levers with parking brake Avid G2 180-mm rotor disc	Avid BB7 mechanical brakes Brake levers with parking brake Gaffer Wave 180mm motor disc	Avid BB7 mechanical brakes Brake levers with parking brake 8 Gear Shimano Alivio Gear Shifter
Wheel	Wheel	Wheel	Wheel	Wheel
12" x 2.40" ultra-grip tyre 12" cast aluminium rims	20" ultra-grip Maxxis 20" black aluminium double-wall rims	18" x 2.40" ultra-grip tyre 18" black aluminium double-wall rims	18" x 2.40" ultra-grip tyre 18" black aluminium double-wall rims	20" Ultra-grip Maxxis 20" black aluminium double-wall rims
Equipment	Equipment	Equipment	Equipment	Equipment
700-lumen front light with position light Dual LED tail lights on the stand frame All-in-one LCD screen with cyclo computer and battery gauge USB charging port Horn Transport handle special for B2	900-lumen front light with position light Dual LED tail lights on the stand frame Display panel and Cycle computr and battery guage, speed regulator and LCD screen USB charging port Removable QR weights	700-lumen front light with position light Dual LED tail lights on the stand frame All-in-one LCD screen with cyclo computer and battery gauge USB charging port Horn Removable QR weights	700-lumen double front light with position light Dual LED tail lights on the stand frame All-in-one LCD screen with cyclo computer and battery gauge USB charging port Horn Removable QR weights	1400-lumen front light with position light Dual LED tail lights on the stand frame Display panel and cycle computer and battery guage, speed regulator and LCD screen USB charging port Horn Removable QR weights
Weight	Weight	Weight	Weight	Weight
12.9kg (without battery)	20.6kg (without battery)	15.9kg (without battery/weights)	15.9kg (without battery/weights)	16.7kg (without battery)

* Autonomy calculated under the following conditions: 75-kg user, ultra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. May vary ±20% depending on the battery.

** Autonomy calculated under the following conditions: 75-kg user, ultra-light wheelchair, flat terrain, no wind and steady speed of 20 km/h. May vary ±20% depending on the battery.

An electric HANDBIKE for every need!

- 1 Discover the full range
- 2 Compare our models
- 3 Choose the one that suits you best!

ALSO AVAILABLE FOR QUADRIPLÉGICS AND HEMIPLEGICS

The versions for **quadriplegics** include the *Gas & Brake* handlebar without levers with which you can accelerate and brake with gentle movements. For **hemiplegics**, the handlebar includes Monolever brakes and the controls can be grouped right or left.

Request a quote or a demonstration today!

Go to our website and request a quote or a demonstration without any commitment. We'll be delighted to organise it for you, as well as to answer any questions you may have.

