

**permobil**  
POWER WHEELCHAIRS



READY SET  
ROLL!

**K300**  
PS Jr.



#### POWERFUL SUSPENSION

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



#### HIGHLY CONFIGURABLE

Allows you to maintain proper positioning as you grow.



#### MULTI-MOTION ARMREST

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

### COLOUR OPTIONS



Volcano Red



Pop Star Pink



Cobalt Blue



Purple Rain



Phantom Grey



Onyx Black

### SPECIFICATIONS

Maximum User Weight	74.8 kg
Maximum Speed	8 kph, 10.4 kph
Range	25.7 km*
Base Width	24" (60.9cm)
Base Length	40" (101.6cm)
Minimum Turning Radius	25" (63.5cm)
Ground Clearance	3" (7.6cm)
Recommended Battery Type	Group 34
Weight Including Batteries	123.3 kg
Drive Electronics	90A R-Net
Maximum Slope	10° Up, 6° Down†
Seat-to-Floor Height	16.5" (elevator) 41.9cm 16.5", 17.25", 18" (seat tube) (41.9, 43.8, 45.7cm)

Optional Seat Elevator	8" (20.3cm)
Tilt Options	0° – 45° (power) 0° – 45° (manual)
Manual Recline	85° – 120°
Legrest Elevation	70° (swing-away) 90° – 150° (center mount elevating)
Backrest Height	16.5" – 22.5" (41.9 - 57.1cm)
Seat Depth	10" – 18" (25.4 - 45.7cm)
Seat Width	11" – 16" (27.9 - 40.6cm)
Armrest Pad Lengths	10", 13" (25.4, 33cm)

\*Actual driving range will vary based on driving and battery conditions.  
†Refer to Owner Manual for safe driving practices.