

**permobil**  
POWER WHEELCHAIRS

ADVENTURE  
AWAITS...

**K450**  
MX

### EXPLORE MORE

Uniquely engineered seating system that brings the world to your fingertips.



### KID-TOUGH TYRES

Large flat tyres improve traction and reduce maintenance.

### ALL-WHEEL SUSPENSION

Suspension system designed for increased balance, improved ride and safety.

High or low, explore more.



### SPECIFICATIONS

Maximum User Weight	56.7 kg
Maximum Speed	8 kph
Range	27.3 km*
Base Width	25" (63.5cm)
Base Length	35" (88.9cm)
Minimum Turning Radius	28" (71cm)
Ground Clearance	2.5" - 6.35cm
Recommended Battery Type	Group 34
Weight Including Batteries	152.8 kg
Drive Electronics	90A R-Net
Maximum Slope	10° Up, 10° Down†
Seat-to-Floor Height	3" - 26" (elevator) [7.6 - 66cm]
Seat Elevator	8" (20.3cm)
Power Tilt	0° - 45°
Manual Recline	85° - 105°
Manual Legrest Elevation	90° - 125°
Backrest Height	14" - 22.5" (35.5 - 57.1cm)
Seat Depth	10" - 18" (25.4 - 45.7cm)
Seat Width	10", 12", 14", 16" (25.4, 30.4, 35.5, 40.6cm)
Armrest Pad Lengths	8", 12.5" (20.3, 31.7cm)

\*Actual driving range will vary based on driving and battery conditions.  
†Refer to Owner Manual for safe driving practices.

### COLOUR OPTIONS



Chili Pepper Red



Pop Star Pink



Extreme Yellow



Titan Blue



Platinum Silver



Onyx Black