

Koala

Ready or not...



permobil




Dynamite performance in a small package means you can safely go wherever the adventure takes you.



Miniflex Power Seating

Unique power seat elevation and tilt designed to provide valuable positioning and exceptional low seat to floor height.



Parental Approved

Optional features such as front bumper and remote stop give any parent that extra piece of mind for their little explorer.

Colour Options



Volcano Red



Popstar Pink



Extreme Yellow



Green Apple



Cobalt Blue



Purple Rain



Platinum Silver



Onyx Black

Additional Hub Colour Options



Chili Pepper Red



Titan Blue



Phantom Gray



Vintage Wine



Crystal Blue



Sublime Green

Specifications

| | |
|----------------------------|------------------------------|
| Maximum User Weight | 35 kg |
| Maximum Speed | 6.5 km/h |
| Range | 24 km* |
| Base Width | 23" |
| Base Length | 32" |
| Minimum Turning Radius | 21" |
| Ground Clearance | 2.5" |
| Recommended Battery Type | M40 |
| Weight Including Batteries | 102.5 kg |
| Drive Electronics | R-Net |
| Maximum Slope | 10° up, 6° down [†] |
| Seat-to-Floor Height | 15" |

| | |
|--------------------------|---------------|
| Optional Seat Elevator | 8" |
| Power Tilt | 0° – 25° |
| Manual Recline | 85° – 120° |
| Manual Legrest Elevation | 90° – 150° |
| Backrest Height | 14", 16", 20" |
| Seat Depth | 10", 12", 14" |
| Seat Width | 10", 12", 14" |
| Armrest Pad Lengths | 8", 10" |

*Actual driving range will vary based on driving and battery conditions.

[†]Refer to Owner Manual for safe driving practices.