

**permobil**  
POWER WHEELCHAIRS

READY OR NOT...




**Koala**  
MINIFLEX



Dynamite performance in a small package means you can safely go wherever the adventure takes you.



**MINIFLEX POWER SEATING**  
 Unique power seat elevation and tilt designed to provide valuable positioning and exceptional low seat to floor height.



**PARENTAL APPROVED**  
 Optional features such as front bumper and remote stop give any parent that extra piece of mind for their little explorer.

**COLOUR OPTIONS**



**SPECIFICATIONS**

Maximum User Weight	36.2 kg
Maximum Speed	6.4 kph
Range	24.1*
Base Width	23" (58.4 cm)
Base Length	32" (81.2 cm)
Minimum Turning Radius	21" (53.3 cm)
Ground Clearance	2.5" (6.3 cm)
Recommended Battery Type	M40
Weight Including Batteries	102.5 kg
Drive Electronics	R-Net
Maximum Slope	10° Up, 6° Down†

Seat-to-Floor Height	15" (38.1 cm)
Optional Seat Elevator	8" (20.3 cm)
Power Tilt	0° – 25°
Manual Recline	85° – 120°
Manual Legrest Elevation	90° – 150°
Backrest Height	14", 16", 20" (35.5, 40.6, 50.8 cm)
Seat Depth	10", 12", 14" (25.4, 30.4, 50.8 cm)
Seat Width	10", 12", 14" (25.4, 30.4, 50.8 cm)
Armrest Pad Lengths	8", 10" (20.3, 25.4 cm)

\*Actual driving range will vary based on driving and battery conditions.  
 †Refer to Owner Manual for safe driving practices.