

**permobil**  
POWER WHEELCHAIRS

HOW SWEET IT IS

**M300**  
PS JR.







Highly adjustable seat that grows as fast as you do.



**CATLINK SUSPENSION**

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



**HIGHLY CONFIGURABLE**

Allows you to maintain proper positioning as you grow.



**MULTI-MOTION ARMREST**

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

**COLOUR OPTIONS**



Volcano Red



Pop Star Pink



Sublime Green



Cobalt Blue



Purple Rain



Platinum Silver



Phantom Grey



Onyx Black

**SPECIFICATIONS**

Maximum User Weight	74.8 kg
Maximum Speed	9.6 kph
Range	20.9 km*
Base Width	24" (60.9cm)
Base Length	36.5" (92.7cm)
Minimum Turning Radius	20" (50.8cm)
Ground Clearance	3" (7.6cm)
Recommended Battery Type	Group 24, Group 34
Weight Including Batteries	141.9 kg G34 / 151 kg G24
Drive Electronics	120A R-Net
Maximum Slope	9° Up, 9° Down†
Seat-to-Floor Height	17" (elevator) 43.1cm 17", 17.75", 18.5" (seat tube) 43.1, 45, 47cm

Optional Seat Elevator	8" (20.3cm)
Tilt Options	0° – 45° (power) 0° – 45° (manual)
Manual Recline	85° – 120°
Legrest Elevation	70° (swing-away) 90° – 150° (center mount elevating)
Backrest Height	16.5" – 22.5" (41.9 - 57.1cm)
Seat Depth	10" – 18" (25.4 - 45.7cm)
Seat Width	11" – 16" (27.9 - 40.6cm)
Armrest Pad Lengths	10", 13" (25.4, 33cm)

\*Actual driving range will vary based on driving and battery conditions.  
†Refer to Owner Manual for safe driving practices.